



DataTrace® Calibration Warranties (U.S. Only)

Mesa Laboratories, Inc. provides a warranty on our NIST traceable calibration services, in addition to our 12 month hardware warranty on new products. To avoid confusion on the Calibration Warranties for DataTrace Tracers, below is a description for the various warranties offered. Note that the factory calibration warranty can vary depending on type of Tracer and whether the subject Tracer is a new build or a Tracer being returned from service.

- All NEW Temperature and Pressure Tracers and the Temperature component of Humidity Tracers are warranted to hold specification for 12 months from date of original shipment.
- All SERVICED Temperature and Pressure Tracers and the Temperature component of Humidity Tracers are warranted to hold specification for 6 months.
- All NEW and SERVICED Humidity Tracers are warranted to hold the Humidity component specification for 3 months.

Field calibration of the humidity sensor is possible on all Humidity Tracers. The Calibration Utility is built into the DataTrace for Windows (DTW) program and the instructions for this function are available in the Operator Manual (Appendix VI) which is on the original DTW CD.

Mesa Laboratories, Inc.
12100 W. 6th Ave.
Lakewood, CO 80228
800-525-1215
303-987-8000